

# HPLC Analysis of Six Miscellaneous Antibacterial Drugs

## Application

### Drug Development

Udo Huber, Adebayo O. Onigbinde

Nalidixic acid, combined with synthetic compounds trimethoprim and sulfamethoxazole, are very effective in treating infections of the urinary tract. Chloramphenicol is a natural compound found in streptomyces, but it can also be easily synthesized. Chloramphenicol is used in treating a wide range of bacterial diseases, for example, typhoid or parrot fever. Furazolidone is used mainly against trichomonades and intestinal bacteria. The analysis of six different antibacterial drugs was obtained with a Zorbax SB-C18 column (Figure 1). Fluorescence spectrophotometry can also be used to detect some of these antibacterial drugs. For example trimethoprim and sulfamethoxazole were analyzed and detected by a fluorescence detector. The fluorescence spectra is shown in Figure 2.




Figure 1. Analysis of six antibacterial drugs using a variable wavelength detector.

## Experimental Conditions


**Equipment:** Agilent 1100 Series HPLC; **UV Detector:** Variable wavelength detector at 0 min, 204 nm, at 4.5 min, 368 nm, at 6 min, 204 nm, at 10 min, 254 nm, at 15 min, 204 nm, standard cell; **Column:** Zorbax SB-C18, 3.5  $\mu$ m, 4.6  $\times$  75 mm (part number 866953-902), Guard cartridges: SB-C18, 5  $\mu$ m, 4.6  $\times$  12.5 mm (part number 820950-920); **Mobile phase:** A = 0.025 M  $\text{KH}_2\text{PO}_4$  in water (pH = 3), B = acetonitrile; **Injection volume:** 5  $\mu$ L; **Temp:** 40  $^\circ\text{C}$ ; **Flow rate:** 1.0 mL/min; **Gradient:** at 0 min 10% B, at 10 min 30% B, at 15 min 60% B; **Column wash:** at 16 min 10% B; **Stop time:** 16 min; **Post time:** 5 min

## Highlights

- The SB-C18 column provides excellent peak shape and selectivity for basic antibacterial drugs.
- The SB-C18 column shows excellent stability at low pH.
- The SB-C18 column shows excellent and rapid resolution of antibacterial drugs at low pH and buffer concentration.
- The HPLC method shows an easy but reliable and precise analysis of the antibacterial drugs.
- The values for limit of detection (LOD), precision of retention time (RT), and area show the good performance of the HPLC analysis.



Agilent Technologies



**Figure 2. Analysis of trimethoprim and sulfamethoxazole using a fluorescence detector (column: Zorbax SB-C18, 5  $\mu$ m, 50  $\times$  2.1 mm)**

Udo Huber is an application chemist based at Agilent Technologies, Waldbronn, Germany.

Adebayo Onigbinde is an applications chemist based at Agilent Technologies, Wilmington, Delaware, USA .

## For More Information

For more information on our products and services, visit our Web site at [www.agilent.com/chem](http://www.agilent.com/chem).

Agilent shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this material.

Information, descriptions, and specifications in this publication are subject to change without notice.

© Agilent Technologies, Inc. 2002

Printed in the USA  
September 10, 2002  
5988-7632EN



**Agilent Technologies**